# CHILL TOWEL



#### DESCRIPTION

THORZT COOLING APPAREL IS DESIGNED TO ASSIST WITH HEAT STRESS PREVENTION BY PROMOTING THE BODY'S THERMOREGULATION.

A cooler body improves the ability to sustain intensive workloads over longer periods of time. The THORZT Chill towel is a great way to cool the skin and contributes to a number of physiological advantages:

- Delay the onset of dehydration
- Less output is directed toward the skin
- A lower core body temperature

The THORZT Chill Towel stays cool for hours with no refrigeration or cold water needed - just 'Wet, Wring and Wave'. It's made from PVA material which can retain moisture for periods far longer than conventional towels, and contains a special plant extract that (once activated) turns moisture cold regardless of it's initial temperature.

### **FEATURES & BENEFITS**

- Promotes Thermoregulation
- Reduces Heat Stress risks
- Comforts and cools the skin
- Looks as cool as it feels
- Machine washable

# PRODUCT DETAILS

Size: 75cm x 32cm Ctn Qty: 20

# CARE INSTRUCTIONS

- It is recommended that the Chill Towel is washed or rinsed before use.
- When finished using, fold the damp Chill Towel and return to packaging do not fold when dry.
- Machine wash with mild detergent.
- Do not put in dryer.





#### CSB

#### **ACTIVATION**

- Soak and Submerge the THORZT Chill Towel in any temperature water.
- 2. Wring out the excess water, leaving it nice and moist.
- 3. Take hold of two corners, and wave it up and down for about 10 seconds.
- 4. Wrap it around your neck or other warm point on the skin for instant cooling relief.

# REACTIVATION

Give it a quick wave every 15-20 minutes to reactivate the cooling process. Re-wet every few hours if needed.

